

DISTRIBUTORS OF CLEANING AND HIGH-LEVEL DISINFECTION MATERIALS

November 2020

Dear Exablate Prostate Customer,

The Exablate Prostate system (Exablate 2100V1, SW Version 8.1) requires Cleaning and High-Level Disinfection after each patient treatment as described in Appendix C of the Operator's Manual (PUB71004685).

This procedure requires materials purchased by the customer, including the following two brand names:

- **Enzol® Enzymatic Detergent**
- **Revital-Ox® RESERT® High Level Disinfectant**

To help locate these materials, Insightec has confirmed that the materials listed above can be ordered and shipped to your country from the distributors listed below.

The Local sales representatives can be found on distributors' websites listed below.

| | |
|--|--|
| ASP (Enzol®): https://www.asp.com/emea/contact-us | STERIS (Revital-Ox®): https://www.steris.com/contact/global-facilities |
|--|--|

In the event you are unable to order the materials listed above, please contact your Insightec representative for assistance.

We look forward to continuing to support you in patient care.

INSIGHTEC

HEADQUARTERS:

5 Nachum Heth St.
Tirat Carmel 39120, Israel
T: +972.4.813.1313

US HEADQUARTERS:

801 Brickell Avenue, Suite 1600
Miami, FL 33131
T: +1.786.534.3849



www.insightec.com